


PulseCath iVAC 2L Post-Market Surveillance Registry Report



	PulseCath iVAC2L Post-Market Surveillance Registry Report	Date: 29-06-2022
		R1058-1

Summary

Retrospective data on the use of PulseCath iVAC2L in the Catheterization Laboratory. Baseline data and clinical outcomes are presented and compared with previous literature on mechanically-assisted high-risk PCI.

1 Abbreviations

AMI:	Acute Myocardial Infarction
CAD:	Coronary Artery Disease
ECMO:	Extra-corporeal Membrane Oxygenation
EF:	Ejection Fraction
IABP:	Intra Aortic Balloon Pump
IQR:	Interquartile Range
LM:	Left Main
MACCE:	Major Adverse Cardiovascular and Cerebral Events
MAP:	Mean Arterial Pressure
PCI:	Percutaneous Coronary Intervention
PMS:	Post-market Surveillance
VA-ECMO:	Veno-arterial Extra-corporeal Membrane Oxygenation

	PulseCath iVAC2L Post-Market Surveillance Registry Report	Date: 29-06-2022
		R1058-1

2 Introduction

Recent technological developments in Interventional Cardiology have enabled PCI in patients with complex coronary artery disease. PulseCath iVAC2L aims to reduce the risk of hemodynamical deterioration during manipulation of the coronary vessels. Post-market surveillance (PMS) evaluates safety and efficacy.

3 Materials and Methods

3.1 Description of iVAC2L. iVAC2L is a percutaneously-inserted left ventricular assist device that ejects arterial blood into the ascending aorta using counterpulsation. Left ventricular blood is aspirated in systole and ejected in the arterial system in diastole, in synchrony with the coronary flow. The aortic pressure waveform during iVAC2L support demonstrates a diastolic “plateau” that partially interrupts the diastolic pressure descent.

3.2 Data Collection. The iVAC2L registry collects clinical data from patients undergoing PulseCath iVAC2L supported high-risk PCI and is derived from published studies, medical records and reports by PulseCath personnel who are on-site during the interventions. The current version of the registry includes cases that occurred in 67 different centers in Europe, South America and Asia.

3.3 Clinical Endpoints. For the analysis of clinical endpoints, only data directly available to the researcher or that has been made available in published peer-reviewed publications are used. Major Adverse Cardiac and Cerebrovascular Events (MACCE) is defined as composite of death (all-causes), acute myocardial infarction, repeat revascularization and Cerebrovascular Events. Intra-procedural Complications was defined as the composite of major bleeding, major vascular complication or respiratory failure. Hemodynamic instability was defined as any situation involving severe hypotension or shock, life-threatening arrhythmias, reported use of vasopressors an/or inotropes, or need to escalate support to another device.

3.4 Data Analysis. Baseline data is presented as means \pm SD or medians \pm IQR as appropriate. Baseline characteristics and clinical outcomes were reported. Clinical outcomes found in the iVAC2L registry were compared with the same endpoints reported in the BCIS-1 and PROTECT II [1, 2] trials using the Fishers’ exact test. Support time and Device Output were compared between the iVAC2L registry and the means / medians from the same trials using Wilcoxon’s one sample test or Student’s t-test for one sample as appropriate. For better comparability with the literature, MACCE was defined as the composite endpoint of death, myocardial infarction, repeat revascularization and stroke. The analysis of clinical outcomes was performed in an intention-to-treat basis, and encompassed all patients in the study. The 2-sided α -error was set to 5% in all comparisons. All calculations were made using the R statistical package version 4.0.2.

4 Results

Data composition. The registry included a total of 174 patients from three previously published studies with iVAC2L (n = 63). Two studies were completed in the Netherlands and a third one in Germany. Another four cases were added from published results. The first described a case of severe acute myocarditis complicated with cardiogenic shock and the second described a case of acute decompensated heart failure. Other two cases described

	PulseCath iVAC2L Post-Market Surveillance Registry Report	Date: 29-06-2022
		R1058-1

successful PCIs for CAD [3-9]. Additional clinical data was collected onsite by PulseCath personnel (n = 108).

Demographics. Demographic characteristics are shown in **Table 3**. Mean age was 69±11 years. Patients tended to have multivessel disease and low ejection fraction (EF). Stenting of the left coronary and its branches was more common than stenting of the right coronary artery. The mean Left Ventricular Ejection Fraction (EF) was 34±14%, and was lower than 40% in 67% of the cases. Significant LM obstruction and three-vessel disease were present in 59% and 54% respectively.

Procedural characteristics. Procedural characteristics are exposed in **Table 4**. The median support time (IQR) was 71 (50-114) minutes, and the average flow produced by iVAC2L was 1.5±0.2L/min. Rotational atherectomy was applied in 11% of the cases. Intraprocedural complications (composite of major bleeding, major vascular complication or respiratory failure) occurred in 10.7% of all the available cases and resulted in the removal of the device in only two cases (1.2%). Hemodynamical instability was reported in 7.7%. Intraprocedural death and MACCE rates accounted for 0.58% and 4.9% of the cases.

Clinical Endpoints. Clinical endpoints are exposed in **Table 4** and put into perspective in **Figure 1**. All-causes mortality rate was 4.6% after 30 days. Acute kidney injury and cerebrovascular events were reported in 6% (each) of 66 individuals who had this information available until discharge. The rate of MACCE after 30 days was 12% out of 66 individuals and 4.9% intraprocedurally.

5 Discussion

This is the largest collection of clinical data on PulseCath iVAC2L to date. In this second year of data collection, this analysis suggest a better safety profile with iVAC2L as compared to other devices based on lower rates of acute myocardial infarction, major bleeding and repeat revascularization.

iVAC2L can also be used in patients receiving Veno-Arterial Extra-corporeal Membrane Oxygenation (VA-ECMO), and as a circulatory backup in severe decompensated heart failure [4, 6, 10, 11]. However, our experience indicates that the most frequent application is during high-risk PCI, being hemodynamically stable at baseline.

The data shows a 12.1% rate of MACCE after 30 days. In contrast, pooled data from two early studies performed in high-risk PCI by Briguori et al [12,13] indicate that the expected MACCE rates in high-risk PCI without MCS may be up to 17%.

With exception of the USPELLA registry, higher rates of MACCE have been reported in two of the main sources of randomized controlled data on MCS in high-risk PCI [1, 2]. The difference, however, did not reach statistical significance.

The Europella registry report, published in 2009 [14], shows a 0.7% rate of hemolysis, higher than the numbers found on the iVAC 2L. In the present analysis, no cases of hemolysis were observed and hence it is reasonable to expect that the real rates would be actually higher than that.

Despite being relatively low, the rate of MACCE with iVAC 2L may have been inflated by the number of cerebrovascular accidents, (CVE), which was proportionally higher in comparison with other studies. The observed rate (6.1%) is due to 4 events that occurred in 3 different sites, three in patients enrolled in the PULSE trial, and one at the Johannes Wesling Hospital (JWH). One patient had a stroke reported in the ICU several days after the intervention, while recovering from a period of refractory ventricular fibrillation. Other two cases consisted in one TIA and one confirmed stroke in two patients treated at the St Thomas Hospital. Both

	PulseCath iVAC2L Post-Market Surveillance Registry Report	Date: 29-06-2022
		R1058-1

individuals had previous history of CVE and one of them was found to have fibromuscular dysplasia of the carotid artery, an uncommon condition that works as an important determinant factor for cerebrovascular events. The fourth subject in the JWH was also reported to have a TIA (thus transitory lasting < 24h) after the intervention.

Importantly, the CVE rate should be regarded in light of the patient profile in the analyzed studies which is characterized by high-risk for cerebrovascular events. Secondly, the studied population in the present report is comparatively smaller than in other studies. A small sample size can compromise data analysis by amplifying sampling error. And thirdly, while it is reasonable to expect that the manufacturer (PulseCath) is to be promptly contacted by operators in face of adverse events that are potentially related to the marketed device, that did not happen in any of the 108 other cases recorded in the PMS registry. If only reported cases are to be considered as occurring events, that would lead to a CVE rate of less than 3% in the entire cohort.

The IMP-IT study [15], which consisted in a retrospectively collected registry of patients receiving any Impella devices for high-risk PCI (n = 117) or cardiogenic shock, provided further information on real-world mortality rates with Impella devices. The time-to-event data released allows to conclude that the mortality rate after 30 days in the high-risk PCI cohort was 9%, which is higher than in this analysis (9% vs 3% with iVAC2L, p = 0.13).

Regarding major bleeding complications, the rate found on the iVAC 2L PMS registry compares favorably with previous data on large-bore catheterization procedures such as transcatheter aortic valve implantation (TAVI), which by default show catheter diameters larger than 7Fr, **Table 4**). The observed rates of vascular complications are most likely reflective of the level of expertise of the operators in the centers involved in the analyzed data which are mostly world-class institutions. The current rate of major vascular complications is higher than in the BCIS-1 and PROTECT II. However, when current rates are put into perspective side by side with relevant studies on MCS and large-bore catheter interventions such as TAVI (**Figure 2**), it becomes evident that several other studies have reported higher rates with different devices, including the IABP (8Fr).

Furthermore, published data suggests that iVAC 2L may be less likely to damage the cardiac structures and to cause hemolysis than Impella. Impella results in damage to the aortic valve in 13% of the cases, while this has never been observed with iVAC 2L. In addition to that, the risk of hemolysis appears to be lower than similar devices. To date, no reports of clinically relevant hemolysis have been received by the manufacturer.

And finally, iVAC 2L apparently succeeds in creating a more stable setting for the coronary interventions to develop. Current data shows a 7.7% rate of hemodynamical instability as opposed to 12% in the BCIS-1 trial, 12.3% in the IABP arm of the PROTECT II study, and 10.2% in the Impella arm of the same study.

iVAC2L induced significant increases in MAP, SBP, DBP and CPO. HR, CO and mPCWP were unaffected. Similarly, an increase in MAP have been described both in the Dutch and in the German studies. The observed increase in CPO is interesting from the clinical perspective because a CPO < 0.6 Watts has been previously related to worse clinical outcomes [16]. Current results consistently indicate that hemolysis levels are lower than in other modalities of mechanical circulatory support [6, 10]. In line with this the present analysis shows no reports of clinically relevant hemolysis.



PulseCath iVAC2L Post-Market Surveillance Registry Report

Date: 29-06-2022

R1058-1




Figure 1. Clinical Endpoints in the iVAC2L PMS Registry and previous data from major studies in Mechanical Circulatory Support.




Figure 2. Clinical Endpoints with PulseCath iVAC2L in the PULSE trial (yellow bars) and in other studies involving MCS and large-bore devices.

	PulseCath iVAC2L Post-Market Surveillance Registry Report	Date: 29-06-2022
		R1058-1

This analysis has the advantage of including a relevant proportion of on-site collected data, reflecting practice in real life. In addition, it also includes data from three prospective studies leading to more standardized definitions and very low levels of missing data. Nevertheless, this comprises only 26% of the registry. Even though there is partial feed-back from operators, this has been mostly addressed by the new Data Collection Forms.

6 Conclusion

In conclusion, this second interim analysis of the iVAC2L PMS registry shows low rates of adverse events and significantly better hemodynamics with iVAC2L. Further research is needed to confirm that in the long run.

7 References

- Perera D, Stables R, Thomas M, et al. Elective intra-aortic balloon counterpulsation during high-risk percutaneous coronary intervention: a randomized controlled trial. *JAMA*. Aug 2010;304(8):867-74. doi:10.1001/jama.2010.1190
- O'Neill WW, Kleiman NS, Moses J, et al. A prospective, randomized clinical trial of hemodynamic support with Impella 2.5 versus intra-aortic balloon pump in patients undergoing high-risk percutaneous coronary intervention: the PROTECT II study. *Circulation*. Oct 2012;126(14):1717-27. doi:10.1161/CIRCULATIONAHA.112.098194
- den Uil CA, Daemen J, Lenzen MJ, et al. Pulsatile iVAC 2L circulatory support in high-risk percutaneous coronary intervention. *EuroIntervention*. Feb 2017;12(14):1689-1696. doi:10.4244/EIJ-D-16-00371
- B Bastos M, van Wiechen M, Van Mieghem N. PulseCath iVAC2L - next generation pulsatile mechanical circulatory support. *Future Cardiology*. 2020;
- A, Schmidt S, Zeise M, Wiemer M, Luani B. High-risk PCI under support of a pulsatile left ventricular assist device - First German experience with the iVAC2L system. *Int J Cardiol*. 12 2019;297:30-35. doi:10.1016/j.ijcard.2019.10.020
- Tschöpe C, Alogna A, Faragli A, et al. Case Report First-in-Man Method Description: Left Ventricular Unloading With iVAC2L During Venous-Arterial Extracorporeal Membrane Oxygenation: From Venous-Arterial Extracorporeal Membrane Oxygenation to EC-MELLA to EC-iVAC®. *Front Cardiovasc Med*. 2020;7:563448. doi:10.3389/fcvm.2020.563448
- Tzikas S, Papadopoulos CH, Evangelidou AP, Vassilikos V. First implantation of the pulsatile left ventricular assist device iVAC2L in a heart failure patient infected with influenza type A. *Hellenic J Cardiol*. May 2020;doi:10.1016/j.hjc.2020.05.002
- Van Mieghem NM, Daemen J, Lenzen MJ, Zandstra R, Malkin O, van Geuns RJ. The PulseCath iVAC 2L left ventricular assist device: conversion to a percutaneous transfemoral approach. *EuroIntervention*. Nov 2015;11(8):835-9. doi:10.4244/EIJV1117A168
- Alekyan BG, Karapetyan NG, Novak AY, Plotnikov GP, Knysh YB, Grigoryan GS. The clinical case of the first high risk percutaneous coronary Intervention under use of the PulseCath iVAC2L left Ventricular assist device in Russia. *Complex Issues of Cardiovascular Diseases* 2020;9(1):103-109 (In Russ). 2020;9(1):103-109.
- Anastasiadis K, Chalvatzoulis O, Antonitsis P, Tossios P, Papakonstantinou C. Left ventricular decompression during peripheral extracorporeal membrane oxygenation support with the use of the novel iVAC pulsatile paracorporeal assist device. *Ann Thorac Surg*. Dec 2011;92(6):2257-9. doi:10.1016/j.athoracsur.2011.05.063
- den Uil CA, Van Mieghem NM, B Bastos M, et al. Primary intra-aortic balloon support versus inotropes for decompensated heart failure and low output: a randomised trial. *EuroIntervention*. 09 2019;15(7):586-593. doi:10.4244/EIJ-D-19-00254
- Briguori C, Sarais C, Pagnotta P, et al. Elective versus provisional intra-aortic balloon pumping in high-risk percutaneous transluminal coronary angioplasty. *Am Heart J*. 2003;145(4):700-707.
- Briguori C, Airoldi F, Chieffo A, et al. Elective versus provisional intraaortic balloon pumping in unprotected left main stenting. *Am Heart J*. 2006;152(3):565-572.
- Sjauw KD, Konorza T, Erbel R, Danna PL, Viecca M, Minden HH, Butter C, Engström T, Hassager C, Machado FP, Pedrazzini G, Wagner DR, Schamberger R, Kerber S, Mathey DG, Schofer J, Engström AE, Henriques JP. Supported high-risk percutaneous coronary intervention with the Impella 2.5 device the Europella registry. *J Am Coll Cardiol*. 2009 Dec 15;54(25):2430-4.
- Chieffo A, Ancona MB, Burzotta F, Pazzanese V, Briguori C, Trani C, Piva T, De Marco F, Di Biasi M, Pagnotta P, Casu G, Giustino G, Montorfano M, Pappalardo F, Tarantini G; Collaborators. Observational multicentre registry of patients treated with IMPella mechanical circulatory support device in Italy: the IMP-IT registry. *EuroIntervention*. 2020 Feb 7;15(15):e1343-e1350.
- Fincke R, Hochman JS, Lowe AM, Menon V, Slater JN, Webb JG, LeJemtel TH, Cotter G; SHOCK Investigators. Cardiac power is the strongest hemodynamic correlate of mortality in cardiogenic shock: a report from the SHOCK trial registry. *J Am Coll Cardiol*. 2004 Jul 21;44(2):340-8.



8 Tables

Country	Number	First Case
UNKNOWN	1	01/01/2019
RUSSIA	1	01/11/2018
TURKEY	1	01/11/2019
JAPAN	1	05/12/2021
CHILE	1	18/08/2021
FRANCE	2	01/01/2019
BELGIUM	2	01/01/2020
SPAIN	2	08/02/2021
BELARUS	2	25/02/2021
FINLAND	4	01/04/2019
GREECE	4	01/11/2019
UAE	4	01/11/2019
SERBIA	4	12/01/2021
HUNGARY	5	01/01/2019
SLOVENIA	5	01/11/2019
INDIA	6	01/01/2020
ITALY	6	01/05/2019
UNITED KINGDOM	7	01/03/2018
POLAND	8	01/04/2019
CROATIA	8	01/12/2020
IRAN	12	01/12/2020
SLOVAKIA	21	01/01/2019
GERMANY	21	01/12/2019
THE NETHERLANDS	46	01/01/2017

Table 1. Number of cases treated with iVAC2L by country and date of the first case.

Table2. Cases by site

Site	Country	Cases
REP. CENTER OF CARDIOLOGY MINSK	BELARUS	2
CHU SART-TILMAN LIEGE	BELGIUM	1
CHU TIVOLI	BELGIUM	1
SOTERO DEL RIO	CHILE	1
CHC OSJEK	CROATIA	1
CHC RIJEKA	CROATIA	1
UHC ZAGREB	CROATIA	3
UNIVERSITY HOSPITAL CENTRE ZAGREB	CROATIA	3
LKS	FINLAND	1
MEILAHTI HOSPITAL	FINLAND	1
LAPLAND CENTRAL HOSPITAL	FINLAND	2
CLINIQUE PASTEUR	FRANCE	2
CHARITE HOSPITAL	GERMANY	1
JW UNIVERSITY HOSPITAL	GERMANY	20
ATTIKON UNIVERSITY HOSPITAL	GREECE	1
ONASSIS MEDICAL CENTER	GREECE	1
IPPOKRATEIO HOSPITAL	GREECE	2
GOKI HOSPITAL	HUNGARY	1
PECS UNIVERSITY	HUNGARY	1
SEMMELWEIS HOSPITAL	HUNGARY	1
SZEGED UNIVEERSITY	HUNGARY	1
UNKNOWN (HUNGARY)	HUNGARY	1
CHRISTIAN MEDICAL COLLEGE VELLORE	INDIA	1
CMC INDIA	INDIA	1
FORTIESCORTS HEART INSTITUTE	INDIA	1
MAX HOSPITAL SAKET	INDIA	1
MEDANTA HOSPITAL	INDIA	2
MASIH HOSPITAL	IRAN	1
MODARES HOSPITAL	IRAN	1
RAJAE HEAR T CENTER	IRAN	1
RASHID HOSPITAL	IRAN	1
VALIASR	IRAN	1
ATIEH HOSPITAL	IRAN	2
RAJEI HEART CENTER	IRAN	3
MARIA CECILIA HOSPITAL	ITALY	1
UNKNOWN (ITALY)	ITALY	1
CLINICA MEDITERRANEA	ITALY	4
UNIVERSITY HOSPITAL KRAKOW	POLAND	1
WARSAW CLINIC	POLAND	1
WOJSKOWY HOSPITAL	POLAND	1
MSWIA WARSAW	POLAND	2
CENTRAL CLINICAL HOSPITAL	POLAND	3
NATIONAL MEDICAL RESEARCH CENTER OF SURGERY	RUSSIA	1
DEDINJE INSTITUTE	SERBIA	1
IKVB DEDINJE	SERBIA	3
KARDIOCENTRUM NITRA	SLOVAKIA	1
PRESOV FACULTY HOSPITAL	SLOVAKIA	1
VUSCH CARDIO CENTER KOSICE	SLOVAKIA	1
MEDISSIMO MEDICAL CENTER	SLOVAKIA	2
SUSCH BANSTRA BYSTRICA	SLOVAKIA	3
MARTIN UNIVERSITY HOSPITAL	SLOVAKIA	6
CINRE HOSPITAL	SLOVAKIA	7
CELJE HOSPITAL	SLOVENIA	1
GH DR FRANC DERGAANC	SLOVENIA	1
GH IZOLA	SLOVENIA	1
MARIBOR	SLOVENIA	1
UMC LJUBLIJANA	SLOVENIA	1

	PulseCath iVAC2L Post-Market Surveillance Registry Report	Date: 29-06-2022
		R1058-1

GOMEZ HOSPITAL	SPAIN	1
HOSPITAL DE SAN JOAN DESPI	SPAIN	1
AMPHIA	THE NETHERLANDS	2
UMC UTRECHT	THE NETHERLANDS	2
AMPHIA BREDA	THE NETHERLANDS	3
ERASMUS MEDICAL CENTER	THE NETHERLANDS	39
MEMORIAL BAHCELIEVER HOSPITAL	TURKEY	1
AL QASSIMI HOSPITAL	UAE	4
ST THOMAS HOSPITAL	UNITED KINGDOM	7
UNKNOWN	UNKNOWN	1

Table 2. Number of cases treated with iVAC2L according to site and country.


	PulseCath iVAC2L Post-Market Surveillance Registry Report	Date: 29-06-2022
		R1058-1

Table 3. Demographic Characteristics of the iVAC 2L dataset.

Variable	Results
Operators	70
Sites	67
Countries	24
Age (years)	69±11 (n = 158)
SYNTAX I score	35±12 (n = 95)
Weight (Kg)	82±17 (n = 118)
Height (cm)	171±15 (n = 112)
Ejection fraction (%)	34±14 (n = 153)
Gender (Male) (%)	79 (n = 142)
EF < 40% (%)	67 (n = 153)
Three-vessel disease (%)	54 (n = 152)
Stented Left Main (%)	59 (n = 153)
Unprotected Left Main (%)	44 (n = 98)
Left Main Equivalent (%)	36 (n = 36)
Stented LAD and branches (%)	70 (n = 132)
Stented LCX and branches (%)	53 (n = 131)
Stented RCA and branches (%)	36 (n = 133)

Operator, site and country are described as counts. Continuous data are exposed as median (IQR) or mean±SD (n = number of observations available). Frequencies are exposed as percentages (n = number of observations available).

Variable	iVAC2L (2022)	BCIS1	USPELLA	PROTECT II (IABP)	PROTECT II (IMPELLA)
Support Time (min)	71 (50-114)	516 (360-1380)***	60 (6-4320)***	504.6±1308.6***	112.2±161.4
Maximum flow (L/min)	1.5 (1.36-1.60)?	---	2.1±0.2***	---	1.9±0.27***
Interruption of Support (%)	1.2	---	---	---	---
Intraprocedural Death (%)	0.6	---	---	---	---
Intraprocedural Complications (%)	10.7	6	---	---	---
Hemodynamic Instability (%)	7.7	12	6.2	12.3	10.2
All-cause Mortality 30-days (%)	4.6	1.3	3.4	6.2	6.9
CVE (%)	6.1	0.7	0.6	1.9	0.5
AMI (%)	3	13*	1.1	13*	17**
Repeat Revascularization (%)	0	1.7	0.6	6.2*	3.2
Major Bleeding (%)	3	3.7	9.7	4.4	12.5*
Major Vascular Complications (%)	6.1	3.3	4	1.4	0.9*
Acute Kidney Injury (%)	6.1	---	2.8	4.2	4.7
MACCE 30-days (%)	12.1	15.6	8	20	14

Table 4. 30-days clinical outcomes of the iVAC2L registry and two major randomized controlled trials in high-risk PCI, the BCIS-1 and PROTECT II trials. MACCE: composite endpoint of death, myocardial infarction, stroke and repeat revascularization after 30 days.

* p < 0.05 versus iVAC2L (2022)
 ** p < 0.01 versus iVAC2L (2022)
 *** p < 0.001 versus iVAC2L (2022)

Variable	Before	During	After	p.value
Heart Rate (bpm)	74.46±14.81	75.22±14.46	73.51±12.64	0.63
CO (L/min)	4.41±1.26	4.84±1.31	4.69±1.21	0.28
SBP (mmHg)	115.24±22.76	122.39±24.07	124.48±24.82	< 0.01
DBP (mmHg)	60.35±15.89	65.25±15.01	64.66±15.69	< 0.05
MAP (mmHg)	78.94±16.04	84.46±15.71	84.79±16.38	< 0.01
mPCWP (mmHg)	16.44±8.66	17.26±9.57	15.56±8.83	0.69
CPO (Watts)	0.72±0.25	0.84±0.3	0.86±0.29	< 0.05

Table 5. Hemodynamic variations with the use of iVAC2L. CO: Cardiac Output. SBP: Systolic Blood Pressure. DBP: Diastolic Blood Pressure. MAP: Mean Arterial Pressure. mPCWP: Mean Pulmonary Wedge Pressure. CPO: Cardiac Power Output. P-values correspond to two-way ANOVA before, during and after support.