# PulseCath iVAC 2L<sup>®</sup>

## Instructions for Use

		Page
en	Instructions for use	2
fr	Mode d'emploi	10
it	Istruzioni per l'uso	19
es	Instrucciones de uso	28
cs	Návod k použití	37
el	Οδηγίες χρήσης	45
ru	Руководство по использованию	54
sk	Pokyny pre používanie	63
sv	Bruksanvisning	71
no	Bruksanvisning	79
pt	Instruções de utilização	87
lt	Naudojimo instrukcija	96
bg	Инструкции за употреба	104
sr	Uputstva za upotrebu	114
hu	Használati útmutató	122

**PulseCath B.V.**

De Corridor 5

3621 ZA BREUKELEN

The Netherlands

Phone: +31 (0) 346 23 0423

E-mail: [info@pulsecath.com](mailto:info@pulsecath.com)

<http://www.pulsecath.com>

## Instructions for use

### PULSECATH iVAC 2L® (REF LV16 and REF LV17)

STERILE: all contents are sterile. **For single use only.** Do not re-sterilize!

Contents:

- LV16 or LV17 catheter with insertion set
- Membrane pump
- Catheter protector
- Extra PTFE Catheter inner tube

Figure 1: Schematic overview of the iVAC 2L® devices. Sizes and diameters are indicated on the product labels.

#### **Device description**

The iVAC 2L devices are designed to provide circulatory support to patients with impaired left ventricular function. There are two types available with catalogue numbers LV17 and LV16.

At heart rates from 60 to 120 beats per minute, the circulatory support provided by the iVAC 2L is 1.0 – 1.5 L/min. The iVAC 2L functions in combination with an Intra Aortic Balloon Pump (IABP) driver.

#### **Intended Purpose**

The iVAC 2L devices are intended for use in patients with impaired left ventricular function who require left ventricular mechanical circulatory support for up to 24 hours.

The iVAC 2L Tip should be positioned in the left ventricular cavity through the femoral artery.

#### **Intended Users**

The iVAC 2L devices are intended to be used by healthcare professionals trained in interventional cardiology procedures.

#### **Intended patient population and medical conditions to be treated**

The range of patient groups that potentially can benefit from the devices is broad, for example patients suffering from Acute Myocardial Infarction (AMI), Cardiogenic Shock (CS), candidates for high-risk Percutaneous Coronary Intervention (HR-PCI) as well as high-risk ablation and mapping.

#### **Indications**

A cardiologist and/or cardiac surgeon will determine if support with the iVAC 2L is the appropriate option. This is based on one or more of the following clinical and/or anatomical criteria.

Clinical criteria:

- Acute coronary syndromes (ACS)
- Chronic kidney disease
- Diabetes mellitus
- Hemodynamic instability
- Left ventricular ejection fraction (LVEF) of <40%
- Previous cardiac surgery

Anatomical criteria:

- Diffuse Coronary Artery Disease (CAD)

- Last patent conduit
- Multivessel disease
- Severe coronary total occlusion
- Severely calcified lesions needing rotational atherectomy
- Unprotected left main coronary disease involving bifurcation

### **Contraindications**

1. Aortic disease: ascending aortic aneurism, severe aortic wall calcifications
2. Aortic valve disease: aortic valve stenosis, aortic valve insufficiency
3. Aortic valve prosthesis
4. Femoral artery stenosis
5. Aneurism of the aorta
6. Thrombus in left ventricle
7. Ejection fraction lower than 10%
8. Right ventricular failure

### **Warnings**

- Read these instructions carefully before use.
- These instructions describe the use of the iVAC 2L in combination with IABP drivers (Datascope 98XT, CS100, CS300 and Cardiosave Hybrid IABP driver, Arrow ACAT 1 and AutoCAT2Wave and AC3). This manual does not replace the IABP driver's manual.
- Be sure that all relevant personnel are adequately trained in using the iVAC 2L and the IABP driver. Personnel involved in using the iVAC 2L require specific training on the procedure steps and the use of the IABP driver with the iVAC 2L. In addition, a supervised procedure guided by a PulseCath representative is part of the training.
- Do not reintroduce the insertion set or PTFE Catheter Inner Tube once removed from the catheter.
- Do not leave the device dormant for extended periods of time in order to prevent formation of thrombi.
- In case of problems, please contact the manufacturer.
- The iVAC 2L is for single use only. Do not reuse or re-sterilize! Reuse or re-sterilization will compromise the mechanical properties of the device which may lead to device failure. As a result, the patient might be injured or die. Reuse or re-sterilization will also create a risk of contamination of the device. Contamination may cause patient infection, illness or death of the patient.
- All components are sterilized by ethylene oxide or gamma radiation.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.
- Patients should be kept sedated during use.
- Assure the femoral artery diameter is sufficiently large prior to inserting the iVAC 2L. For femoral arteries with a diameter of  $\geq 6.8$  mm the LV17 type is recommended, for femoral arteries of  $\geq 6.0$  mm, the LV16 is recommended.
- Do not use alcohol containing detergents to clean or disinfect the device prior to removal, as these might damage the device which could lead to leakage of parts of the device.
- Do not use tie-wraps to fasten the connection between the catheter and the membrane pump, as tie-wraps can cause leakage of the connection.
- Beware of needle sticks into the catheter: a puncture of the catheter will immediately cause aspiration of air into the device, resulting in the air being ejected in the aorta of the patient. The IABP driver should be switched off immediately when a problem is suspected.
- For Datascope IABP driver, "IAB" in help screens should be read as "iVAC 2L".
- Make sure the membrane pump driveline remains connected to the IABP driver. Also make sure the membrane pump driveline is not kinked or compressed. A loose connection or kink will stop the pumping action of the iVAC 2L.
- It is recommended to flush both ports of the LV16 / LV17 catheter with insertion set every 5 minutes with heparinized saline after insertion, as long as the insertion set has not been taken out. Make sure the insertion set is fully de-aired before flushing.

Accordingly, PulseCath will not be responsible for any direct, incidental or consequential damages or expenses resulting from use by untrained personnel or reuse of the product.

## **Precautions**

- Store in a dry, dark and cool (0 °C – 25 °C) place.
- Do not use open or damaged packages.
- Use prior to the 'Use by' date
- Use a pigtail catheter and corresponding guide wire to introduce the iVAC 2L.
- The product has been tested and qualified with accessories (see necessary equipment and disposables). The use of any other accessory could result in complications and/or malfunction of the iVAC 2L.
- Do not leave the patient unattended during use of the iVAC 2L.
- To prevent thrombosis and product malfunction the patient coagulation has to be suppressed continuously. An ACT of 200 seconds minimum is recommended. Anticoagulation should be monitored regularly.
- Arterial pressure and ECG activity should be monitored continuously during iVAC 2L use.
- Oxygen saturation of the leg should be monitored.
- The iVAC 2L performance is optimal at 1:1 frequency. At 1:2 frequency the performance decreases with 50% and at 1:4 frequency with 75%.
- Keep the mean blood pressure above 60 mmHg in order to obtain optimal circulatory support.
- At any sign of improper positioning of the iVAC 2L (i.e. bad filling, sharp peaks in pressure graph) verify correct position of the catheter by transesophageal echo (TEE) or X-ray (fluoroscopy).
- If severe difficulty or strong resistance is met during any stage of the procedure, discontinue the procedure and determine the cause before proceeding.
- When using an Arrow ACAT1, AutoCAT2Wave or AC3 IABP driver the option to permanently disable the Gas Alarms needs to be enabled by the Distributor of these devices.
- The sheathless introduction must be limited to patients with favorable anatomy. This includes patients without significant atherosclerosis or calcification at the access site and absence of significant disease or tortuosity in the iliac and femoral arteries.

## **Limitations**

- The iVAC 2L is permitted for use up to 24 hours. There is lack of clinical data to support the application of the iVAC 2L beyond this time point.
- The iVAC 2L functions are sub-optimal at heart rates lower than 60 bpm or higher than 120 bpm.
- Use of the internal triggering mode will decrease the functioning of the iVAC 2L.
- The iVAC 2L is not suitable for mobile use.

## **Complications**

Invasive procedures should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during the procedure.

Many of the risks related to an iVAC 2L supported procedure are the same as those with the heart / coronary treatment being completed and the placement of any pump used to help the heart.

Possible complications that have occurred (rate of occurrence) and may occur (not reported) include, but are not limited to the following:

- Acute kidney Injury (<3%)
- Acute Myocardial Infarction (<2%)
- Cardiac arrest requiring CPR (<3%)
- Cerebro-vascular Accident (CVA) (<3%)
- Death (<6%)
- Hemolysis (not reported)
- Induction of mitral valve incompetence (not reported)
- Infection (not reported)
- Injury to the aortic valve (not reported)
- Ischemia of the leg due to obstruction of the femoral artery (<2%)
- Major Bleeding (<2%)
- Mechanical trauma to the heart and inner structures (not reported)
- Perforation of the access site (<2%)
- Repeat revascularization (not reported)
- Severe hypotension (<12%)

- Thrombosis at the insertion site (assuming angiographic finding) (<2%)
- Thrombosis when the device is stopped for a longer period of time (not reported).

Contact PulseCath and the competent authority of the Member State in which the user and/or patient is established in case any serious incident has occurred in relation to the device. Refer to the front page for PulseCath's contact information.

### **Necessary equipment and disposables**

- IABP driver.
- Introducer Sheath, with minimum inside diameter of 17 Fr for the LV16 catheter and 18 Fr minimum for the LV17 catheter.
- Guide wire: 0.035" or 0.038", length 260 cm (Super Stiff) and corresponding needle.
- Heparinized saline (2500 IU heparin in 500 mL saline).
- When using an Arrow ACAT1, AutoCat2Wave or AC3 IABP driver: Arrow Pump adapter for 50 cc IABs (Datascope, ref. 0684-00-0501-02).

### **Preparation of the patient**

Additionally to normal clinical proceedings it should be assured that:

- The patient is sufficiently heparinized, an ACT of 200 seconds minimum is recommended during use of the iVAC 2L.
- Arterial pressure and ECG are monitored continuously.
- Arterial pressure and ECG signals are connected to the IABP driver.
- The diameter of the entry artery is measured to determine whether it is sufficiently larger than the iVAC 2L catheter diameter.
- Oxygen saturation of the leg is monitored as a control on peripheral perfusion.

### **Preparation of the iVAC 2L devices**

- Open the pouches in such a way that the devices remain sterile.
- The PTFE Catheter Inner Tube is supplied with a protective piece of tubing on the tip. Remove the piece of tubing before proceeding.
- Exchange the Catheter inner tube with the extra inner tube available in the box by pulling out the existing inner tube and inserting the Extra PTFE Catheter Inner Tube carefully until it protrudes the Catheter tip.
- Flush the inner lumen of the catheter and the inner lumen of the insertion set with heparinized saline through the sideline extending from the plug at the proximal end of the catheter and the sideline of the hemostasis valve at the proximal end of the insertion set.
- Make sure that all stopcocks are in the open position to let escape the air during insertion.
- De-air the membrane pump by filling it with heparinized saline and remove all air-bubbles by shaking and by ticking against its housing.

### **Preparation of the IABP driver**

- Switch on the IABP driver and open the helium gas bottle.
- For Datascope IABP drivers: make sure an Adult safety disc is used.
- For Arrow IABP drivers: make sure an Arrow Pump adapter for 50cc IABs is present.

As the iVAC 2L has a different resistance than an IAB, additional actions, depending on the type of IABP driver should be taken:

### **Datascope 98XT IABP driver:**

1. Set Trigger Select on ECG or AP.
2. Disable the Augmentation Alarm.
3. When the iVAC 2L is connected: press IAB Fill button (press 2 seconds).
4. Let IABP driver fill the membrane pump.
5. Set IAB Fill mode at Manual Fill. The Slow Gas Loss Alarm is now disabled.
6. Disable the R-trac option.
7. Start the IABP driver at 1:2 frequency with maximum augmentation and observe the movements of the membrane of the Membrane pump.

8. When the membrane is moving smoothly and the timing is correct, set the frequency at 1:1 for maximum performance.
9. In case it is necessary to use internal triggering: remove the ECG cable from the IABP driver. Replace the ECG cable when returning to the ECG or AP triggering mode.

**Datascope CS100 and CS300 IABP driver:**

1. Set Operation mode on Semi-Automatic.
2. Set Trigger Source at ECG or AP.
3. Disable the Augmentation Alarm.
4. Disable the R-trac option.
5. When the iVAC 2L is connected: press the Start button.
6. Let the IABP driver fill the membrane pump.
7. Start the IABP driver at 1:2 frequency with maximum augmentation and observe the movements of the membrane of the membrane pump.
8. When the membrane is moving smoothly and the timing is correct, set the frequency at 1:1 for maximum performance.
9. In case it is necessary to use internal triggering: remove the ECG cable from the IABP driver. Replace the ECG cable when returning to the ECG or AP triggering mode.
10. In case of repeating alarms: set Operation mode on Semi-Automatic; Press IAB Fill button (press 2 seconds); Let the IABP driver fill the Membrane pump.
11. Then Set the IABP Fill mode at Manual Fill: Press the Pump options button, scroll through menu and set the Fill mode to manual. The Slow Gas Loss Alarm is now disabled; Press the Start button.

**Getinge/Datascope/Maquet Cardiosave Hybrid IABP driver:**

1. Set OPERATION MODE on SEMI AUTO.
2. Set TRIGGER on ECG or PRESSURE.
3. Disable the AUGMENTATION ALARM (AUG. ALARM).
4. Disable the R-TRAC option (Select PREFERENCES -> PUMP OPTIONS -> R-trac).
5. When the iVAC 2L is connected: press the Start button.
6. Let the IABP driver fill the membrane pump.
7. Start the IABP driver at 1:2 frequency with maximum augmentation and observe the movements of the membrane pump.
8. When the membrane is moving smoothly and the timing is correct, set the frequency at 1:1 for maximum performance.

To set the IABP driver in INTERNAL mode: select INTERNAL in TRIGGER and remove the ECG cable from the IABP driver. Reconnect the ECG cable when returning to the ECG or PRESSURE triggering mode.

In case of repeating alarms: Select PREFERENCES -> PUMP OPTIONS -> CATHETER ALARMS

**Arrow ACAT1 / AutoCAT2Wave / Arrow AC3 IABP driver:**

1. Set the driver in Operator Mode (for AutoCat2Wave driver only).
2. Set the Trigger mode at Pattern, Peak or AP.
3. Set Alarms Permanent Off: press the Alarms button, press Permanent Off and confirm. The Gas Loss alarms are now disabled. The option of permanent disabling the Alarms needs to be enabled by the Arrow distributor.
4. Connect Arrow Pump adapter for 50cc IABs to the membrane pump's driveline. See the Arrow manual for further information on using the Arrow Pump adapter.
5. When the iVAC 2L is connected: start the IABP driver at an 1:2 frequency and observe the movements of the membrane of the Membrane pump.
6. When the membrane is moving smoothly and the timing is correct, set the frequency at 1:1.
7. To restart pumping after a stop: press the Reset button twice and then press the Start button.

### **Introduction procedure**

Introduction of the iVAC 2L into the femoral artery can be performed via an introducer sheath, and using a guide wire to guide the catheter through the aortic valve into the left ventricle. For sheathless introduction of the iVAC 2L into the femoral artery, the same steps apply without the initial introduction of the sheath. Consider positioning of a safety wire to secure a bailout pathway according to accepted clinical practice.

To determine the correct position of the iVAC 2L (tip in left ventricle and valve in aorta) it is strongly recommended to use X-ray (fluoroscopy) or TEE imaging.

The arrow on the Connector of the iVAC 2L catheter indicates the position of the opening of the valve of the iVAC 2L.

1. Prepare the entrance site in the femoral artery according to the hospital standard procedure.
2. Introduce the introducer sheath in the entrance site according to the instruction for use of the sheath.
3. Forward the guide wire into the femoral artery until the tip is in the left ventricle (TEE or X-ray control).
4. Insert the proximal end of the guide wire into the tip of the iVAC 2L catheter, so that the proximal end of the guide wire extends out of the proximal end of the iVAC 2L catheter.
5. Close the haemostasis valve of the iVAC 2L catheter by rotating the cap until there is no bleeding. Leave freedom to move the iVAC 2L catheter over the guide wire.
6. Guide the iVAC 2L catheter over the guide wire into the insertion sheath in the femoral artery.
7. Forward the catheter and de-air the catheter using the sideline extending from the plug at the proximal end.
8. Carefully guide the tip of the iVAC 2L catheter into the ventricle. The tip of the PTFE Catheter Inner Tube protrudes approximately 3 cm from the iVAC 2L catheter's tip. Consider retracting the PTFE Catheter Inner Tube back into the catheter's metal tip before crossing the aortic valve and positioning of the catheter's tip inside the left ventricle.
9. Inspect the position of the iVAC 2L catheter tip by TEE or X-ray.
10. When the position is correct, remove the guide wire.
11. Pull back the iVAC 2L Catheter inner tube until the tip of the tube is in the plug at the proximal end of the iVAC 2L catheter.
12. Place a Tube Clamp in the middle of the Connector of the iVAC 2L catheter.
13. Remove the plug with the iVAC 2L Catheter inner tube.
14. Fill the iVAC 2L catheter Connector, and the membrane pump, completely with heparinized saline and connect them while continuously adding heparinized saline to prevent air-entrapment.
15. Remove the Tube Clamp and make sure (visually) that there are no air bubbles in the membrane pump.
16. In case of air bubbles in the membrane pump, go back to step 13 and disconnect the membrane pump. De-air the membrane as mentioned at "preparation of the iVAC 2L" chapter and return to step 15.
17. Connect the driveline of the membrane pump with the IABP driver.
18. Start pumping at 1:2 frequency with maximal augmentation, see "Preparation of the IABP driver".
19. Adjust timing, see "Operating the IABP driver".
20. When the membrane of the Membrane pump is moving smoothly and the timing is correct: set the frequency at 1:1.

### **Operating the IABP driver**

The control of the IABP driver with the iVAC 2L is similar to the control of the IAB. ECG signal and Aortic Pressure (AP) signal should be connected to the IABP driver. See IABP driver manual for further instructions on handling the driver.

The IABP driver should be set at ECG triggering or at AP triggering. The Internal triggering mode should only be used in case of a very irregular heart rate or at very high heart beat rates.

The iVAC 2L should eject during diastole and aspire during systole. The timing should be set in such a way that the inflation may come as soon as the slope of AP curve is decreasing, just prior to the dicrotic notch. The deflation should be set prior to the systole of the heart. Correct timing should be determined by the shape of the AP curvature.

Smooth movement of the membrane pump indicates unrestricted blood inflow. Disturbed membrane movement or membrane pump vibration during aspiration shows restricted blood inflow. This can be

solved by gently pulling the iVAC 2L catheter backwards. Always check if the tip of the catheter is still in the ventricle.

The “ticking” of the iVAC 2L valve is audible. One “tick” per beat indicates unrestricted functioning of the iVAC 2L. A series of ticks per beat indicates restricted movements. This can be solved by rotating the iVAC 2L catheter.

The shape of the “Balloon Pressure Waveform” gives an indication of correct functioning (Fig 2). A round shape of the pressure peaks indicate correct functioning. Sharp peaks indicate an inflow obstruction, an incorrect position, or bad timing.


Figure 2: Balloon pressure waveform. Left: correct (round) peaks; Right: incorrect (sharp) peaks.

### **Explantation procedure**

During explantation of the iVAC 2L the patient should be sedated. All necessary steps should be undertaken to prevent infection.

When using disinfectants, do not use alcohol-based disinfectants as these can damage the device.

1. Wean the patient by setting the frequency to 1:2 and to 1:4 for a period of time depending on the condition of the patient.
2. Stop the IABP driver.
3. Place a Tube Clamp on the Connector and disconnect the membrane pump.
4. Pull the iVAC 2L backwards. Pull gently, do not put any force on the catheter. Make sure the catheter was removed completely.
5. Close the artery and close the wound.

When the iVAC 2L catheter is introduced without a sheath, it is recommended to reinsert a guidewire through the iVAC 2L catheter to secure the femoral access to be used in case of complications.

**NOTE:** Dispose contaminated products and packaging materials utilizing standard hospital procedures and universal precautions for bio-hazardous waste.

### **Summary of Safety and Clinical Performance**

The Summary of Safety and Clinical Performance (SSCP) of the iVAC 2L devices is available on the website of the manufacturer (<http://www.pulsecath.com>).

-  = Refer to accompanying instructions for use.
-  = Use-by date
-  = Catalogue number
-  = Do not re-use
-  = Lot number
-  = Unique Device Identification
-  = Do not re-sterilize
-  = Sterile product
-  = Sterilized by ethylene oxide
-  = Sterilized by gamma radiation
-  = Single sterile barrier system
-  = Single sterile barrier system with protective packaging outside
-  = Do not use if the product sterile barrier system or its packaging is compromised
-  = Keep away from sunlight
-  = Keep dry
-  = 25°C upper limit storage temperature, 0°C lower limit storage temperature. Transient distribution temperatures may exceed these limits.
-  = Medical Device
-  = Date of manufacture

Manufacturer:

**PulseCath B.V.**

De Corridor 5  
3621 ZA BREUKELEN  
The Netherlands

Phone: +31 (0) 346 23 0423

E-mail: [info@pulsecath.com](mailto:info@pulsecath.com)

<http://www.pulsecath.com>

**DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY**

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the PulseCath product(s) described in this publication. Under no circumstances shall PulseCath be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind PulseCath to any representation or warranty except as specifically set forth herein. Descriptions or specifications in PulseCath printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

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