



Components		
#	Description	Materials
A	Tip	Stainless steel
B	Valve	Stainless steel
C	Flexible catheter segment	Nitinol wire encapsulated in Polyurethane
D	Soft connector	Polyurethane
E	Membrane Pump	MABS, Polycarbonate, Polyurethane
F	Membrane Pump Driveline with Male Luer Slip	PVC (driveline tube), Acrylic compound (Male Luer Slip)

Dimensions		
#	Description	Measurement (mm)
g	Maximum protrusion length of PTFE stylet (dilator) outside catheter tip	30
h	Length between tip and valve	60
i	Catheter diameter	LV16: 5.33 LV17 5.66
j	Catheter length excluding soft connector	920
k	Total length catheter	1000